



COMMUNICATIONS FOR TODAY ... AND TOMORROW.

Tek-CARE500 TekTip: Software Update

Tek-CARE500 vers. 1.3.0.462 for upgrading vers. 1.1.0.349 & above

In a continuing effort to address the issues brought to our attention by customer feedback, a software update has been released for the Tek-CARE500 Wireless Emergency Call System. This revision allows greater flexibility in setting supervision times vs. timeouts, to maintain an optimum level of traffic and to reduce the frequency of transmitters erroneously reporting as missing. The patch also improves database maintenance.

A script file is also now available to edit supervision timeouts universally, thus eliminating the need to change each device individually. To obtain this script file, the software update, or to inquire about updating a Tek-CARE500 system older than version 1.1.0.349, contact Technical Support at 800.332.6608.

(2nd Quarter 2006)

www.tektone.com

324 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net
Phone: 828.524.9967 • Fax: 828.524.9968 • Sales: option 2 • Tech Support: option 3

