

Features

- UL® 1069 Listed
- Supervised
- Easy bulb replacement
- Excellent visibility
- Lamp supervision
- Interchangeable lens caps
- Simplified multiplex wiring
- Plug-in connectors

Description

TekTone's LI386C Dome/Zone Light with built-in multipurpose station provides visual indication of calls originating from bath, emergency and code call switches when used as a dome light. The LI386C also indicates calls originating from patient, staff and duty stations when used as a zone light. The LI386C also indicates all calls from the zones to which those stations are assigned and provides communication to and from the central equipment, including the nurse master station and monitor. It is used with remotely-mounted devices such as switches, contact closures, or other signal-originating equipment. No controls are provided.

Listings



UL® 1069 listed; UL® File Number E71599.



LI386C Addressable Dome/Zone Light

Specifications

Dimensions:	Height:	4.5625"	(115 mm)
	Width:	4.5625"	(115 mm)
	Depth:	3.3"	(84 mm)
Bulbs:	#1829 28v. 70mA		
Power:	24 VDC, 51 mA		
Connections:	Ramp lock plug-in type		
Finish:	Base is flame retardant (UL® 94V-0) plastic. Amber, red and green lenses are flame-retardant translucent polycarbonate plastic. White lens is flame retardant translucent polystyrene plastic.		
Housing:	Standard dual gang box and ring.		

Required Components

Tek-CARE®300II Nurse Call System