

TekTone® Integrated Wired and Wireless Call Systems w/ Domes

IL1091 • Rev. 0 • 02/2022

Tek-CARE120/160
Nurse Call System (Wired)




NC415AV
Master Station



NC120/NC160
Central Equipment (CE) Modules



NC554/8
Network Switch



Tek-CARE500
Wireless Call System



NC511UL
Repeater/Locator



NC511UL
Repeater/Locator



NC510UL
Receiver



SF521UL/SF523UL
Wireless Stations



SF503UL/SF504UL
Pendant Transmitters

Facility
LAN



NC475DESK
Tek-CARE® Appliance Server

Wiring

See individual layout sheets/manuals
for complete wiring information.

Tek-CARE120	IL1004/IL992
Tek-CARE160	IL1025/IL1024
Tek-CARE500	IL814/IL881