
Description

The IH431S is a surface-mounting backbox for IR43xG3 stations. Stations snap-fit to the box so no screws are needed to mount stations to the backbox. The IH431S is constructed of durable ABS plastic, and mounting holes align with 2- and 3-gang electrical boxes.

Required Components

Tek-CARE®400 GEN3 Nurse Call System
IR43xG3-series Patient Stations

Specifications

Dimensions: 5"H × 7.56"W × 2.6"D
(127 mm × 192 mm × 66 mm)

Construction: ABS plastic

Mounting: Screw holes align with 2- and 3-gang electrical boxes



IH431S Surface Backbox
