

Tek-CARE500 vs. Tek-CARE570 System Comparison

Our newest Wireless Nurse Call System, The Tek-CARE570, was specifically designed for smaller single building facilities. However, our Tek-CARE500 Wireless Nurse Call System was designed for facilities in a campus environment. Several other factors determine which system is a best fit. Find out which Wireless Nurse Call System is right for your facility with our comparison chart below.

	Tek-CARE500 	Tek-CARE570 
Intended Application	Campus Environment	Single Building
UL Listing	UL2560	UL2560
Water Resistant Wearable Pendants	✓	✓
Water Resistant Wall-Mounted Emergency Stations	✓	✓
Fully Supervised	✓	✓
Resident Information Page	✓	✓
Resident Check-in	✓	✓
Door/Window Monitors	✓	✓
Smoke Detectors	✓	
Motion Detectors - Check-in	✓	
Maximum Number of Transmitters	3000	900
Tek-ALERT Integration	✓	✓
Reporting	✓	✓
Mobile App	✓	✓
Dome Lights *	✓	✓
Master Map View	✓	✓
Enhanced Location Capability **	✓	pending

*Requires use of Tek-CARE120

**TEP only