

### Features

- Tactile push-button call cord
  - Call button makes an audible and tactile click for assurance of activation
- UL® 1069 Listed
- Oxygen-safe
- Impact-resistant and waterproof
- Small, lightweight pendant fits comfortably in hand
- Bed clip included
- 8-pin DIN plug

### Description

The TekTone SF401B call cord is for use with the Tek-CARE3000, Tek-CARE400 P5/P5+/GEN3, and Tek-CARE300III Nurse Call Systems, where 8-pin DIN jacks are used. The SF401B is an oxygen-safe call cord that is certified safe for use by all patients receiving oxygen therapy.

This waterproof and impact-resistant call cord has a completely sealed push button and a polyvinyl-covered cord.

### Required Components

Tek-CARE®3000 Nurse Call System with IR331 or IR332 Patient Stations

Tek-CARE400® GEN3 Nurse Call System  
IR430G3, IR431G3/IR432G3, IR434G3 patient stations

Tek-CARE400® Nurse Call System  
IR411, IR412, IR412EPB patient stations

Tek-CARE400® P5/P5+ Nurse Call System  
IR411P5/IR421P5, IR412P5/IR422P5 patient stations

Tek-CARE®300III Nurse Call System  
IR319B, IR320B patient stations



SF401B Call Cord

### Specifications

Dimensions:	SF401B:	8' long
	SF401B/10:	10' long
	SF401B/20:	20' long
Connection:	8-pin DIN plug	
Construction:	SPST1-A switch Integral brass, fine silver contacts 1 AMP, 15 watt maximum	

### Listings



UL®1069 Listed and cUL® Listed to CSA  
C22.2 No.205; UL® File Number E165723.