

### Features

- Routine audible-visual signaling
- Emergency audible-visual signaling
- Voice communication and monitoring
- Push button operation
- Nurse follow and group call features
- Compact modular design for wall mounting
- Integrated circuit technology

### Description

The basic Tek-CARE®NC200 series Audio/Visual Nurse Call System consists of patient calling stations; emergency call stations; a master intercom/annunciator station including amplifier, dome lights, and zone lights; and reliable, solid-state control and power supply equipment. The flexible, modular design permits simple, customized installation.

Emergency signals are distinctively identified by a rapid pulsating of all visual and audible signals associated with the calling station. If routine and emergency calls are resent simultaneously, emergency signaling is given precedence. Emergency audible signals can only be canceled at the point of origin.

Sensitive two-way voice communication, Nurse Follow, Tone-Off, Tone Call, and Group Call functions are push button selected.

Optional equipment is available to provide for staff stations, duty stations, remote visual annunciator, master handset, and a third priority level annunciation.

### Listings



UL® 1069 listed; UL® File Number E71599.



NC200N Master Panel with 2 PM200/20N Button Panels, and TA200N Handset

### Required Components

NC200N Master Station Panel  
PM200/20N Add-on Panel with 20-station selector/light buttons, test button, and group call buttons  
IH151NK Power and Control Kit

### Optional Components

TA200N Handset  
PK124K Battery Backup Kit  
PK153 Third Priority Control Unit