

Features

- Simple, push-button calling using plug-in call cords
- Combined call-placed indicator and cancel button
- Call automatically registers if either call cord is accidentally removed
- Flame retardant plastic faceplate
- Corridor zone light capability

Description

The SF102 Dual Bed Room Station includes two jacks allowing the attachment of two SF301-series Call Cords. A call is registered at the nurse station, corridor light and corridor zone light (if installed) when either cord button is depressed or if either call cord is accidentally removed from the jack on the SF102.

The station includes a lighted call-placed switch that lights when a call is placed. The call-placed switch also provides the system cancel function. All of the electronics are mounted on an attractive, flame retardant plastic faceplate.

Listings



UL 1069 Listed; UL® File Number E71599.

Specifications

- Dimensions: 4.5"H × 2.75"W × 1.75"D
(114 mm × 69 mm × 44 mm)
- Connections: Plug-in connectors
- Construction: Flame retardant ABS plastic faceplate.
Switch and electronic components are circuit board mounted
- Housing: Single-gang ring, or single-gang ring with dual-gang box



SF102 Dual Bed Room Station

Required Components

Tek-CARE NC110 Nurse Call System
SF301-series Call Cords