



## Tek-CARE®120 & Tek-CARE®160 Hardware Matrix

### Control Modules

Part	Module Description	Function	Spec Sheet
NC120	Gen2 CE Module	Controls master and station communication on the Tek-CARE120	IL994
NC160	Gen2 CE Module	Controls master and station communication on the Tek-CARE160	IL1014
NC475	Tek-CARE Appliance Server	Provides control for NC404TS master station and serves as platform for all software add-ons (Tek-CARE Staff App, Reporting, etc.), wireless telephone interface and third-party integrations	IL1011

### Master Stations

Part	Display	Non-Touchscreen Controls	Spec Sheet
NC404TS	22" color touchscreen display	Calls, Faults, Talk, and Reset button on NC404 side cradle handset	IL1034
NC415AV	5.5" color touchscreen display	Talk, Reset button, optional TA415 handset	IL1058

## Stations

Part	Station Type	Series	Front Inputs	Spec Sheet
SF121	Single Patient station	120,160	1/4" jack, reset button	IL998
SF122	Dual Patient station	120,160	two 1/4" jacks, reset button	IL998
IR160	Audio station/remote speaker	160	N/A	IL1013
SF123	Customizable 2-Button Pull Cord	120,160	Customizable push button, reset button, optional pull cord	IL999
SF125	Duty station, 4 annunciation LEDs	120,160	N/A	IL1000

## Room Controllers

PM120	Room Controller	120,160	N/A	IL997
LI122(BL)	Room Controller/Dome Light	120,160	N/A	IL996