

### Features

- Non-powered pneumatically operated switch
- Large 3.5" diameter pad for easy activation
- Oxygen safe
- Latex free
- Includes security clip

### Description

The SF401GP Geriatric Pad Call Cord is a round, flattened soft vinyl pad with gripping feet to prevent it from slipping off the patient's bed. The SF401GP works well for residents who have limited hand dexterity. It has a mechanical switch that is activated by slight pressure from the hand, arm or body, or can even be positioned under the resident's chin.

The SF401GP has an 8-foot cord and an 8-pin DIN plug. Also available is the SF401GP/10 with a 10-foot cord.

### Required Components

Tek-CARE®400 P5/P5+/GEN3 Nurse Call Systems  
IR4xxP5-series or IR43xG3-series stations  
SF401EX Strain Relief Cable

Tek-CARE®300III Nurse Call Systems  
IR319- or IR320-series Patient Stations  
SF4-series Call Cords or Pillow Speakers

Tek-CARE®3000 Nurse Call System  
IR331 or IR332 patient stations  
SF401EX Strain Relief Cable



SF401GP Geriatric Pad Call Cord

### Specifications

Cord length:	8' (2.44 m)
Connections:	8-pin DIN plugs
Switch:	Momentary contact
Electrical rating:	24 volts AC or DC at 0.33 amps

### Listings



UL® 1069 Listed and cUL® Listed to CSA  
C22.2 No. 205; UL® File Number E71599.