

Features

- Large red nurse call button
- Three AUX call buttons for Water, Pain, and Toilet
- Four environmental control buttons
- 8-foot cord with DIN plug and security clip
- Provides pillow speaker intercom with built-in microphone, environmental control, and television control
- UL® 1069, and UL® 2560 Listed
- Oxygen safe
- Fully sealed for easy cleaning
- Flame-retardant (UL® 94V-0) high-impact plastic

Description

The SF431PSG3 Pillow Speaker is used with the Tek-CARE®400 GEN3 IR431G3 and IR432G3 patient stations, and the Tek-CARE®3000 IR331 and IR332 patient stations.

The pillow speaker provides call placement and audio communication to the Tek-CARE System and offers remote control for hospital-grade TVs.

The pillow speaker also supplies two-way, hands-free communication with the master station, and automatic silencing of TV audio while the intercom is in use. Water, Pain, Toilet, two light control buttons (typically ROOM and READ lights) and two additional environmental controls (ex. temperature up/down).

Push buttons are provided to call the nurse, plus TV power, channel selection up/down, Closed Captioning, and mute.

Required Components

Tek-CARE®400 GEN3 Nurse Call System
IR431G3 or IR432G3-series Patient Station

Tek-CARE®3000 Nurse Call System
IR331 or IR332 Patient Station

Optional Components

Low Voltage Controller for light and temperature control
(supplied by others)



SF431PSG3 Pillow Speaker

Listings



UL® 1069, UL® 2560, and cUL® listed to CSA
C22.2 No. 205; UL® File Number E71599.